

**NOTE: Revision to  
Assay Protocol**

---

**Product Manual**

# **OxiSelect™ Comet Assay Kit (3-Well Slides)**


**Catalog Number**

**STA-350**

**15 assays**

**FOR RESEARCH USE ONLY  
Not for use in diagnostic procedures**

---



**CELL BIOLABS, INC.**  
*Creating Solutions for Life Science Research*

## **Introduction**


DNA damage, due to environmental factors and normal metabolic processes inside the cell, occurs at a rate of 1,000 to 1,000,000 molecular lesions per cell per day. While this counts for only a small part of the human genome's approximately 6 billion bases (3 billion base pairs), unrepaired lesions to critical genes can impede a cell's ability to carry out its function and appreciably increase the likelihood of cancer.

The comet assay, or single cell gel electrophoresis assay (SCGE), is a common technique for measurement of DNA damage in individual cells. Under an electrophoretic field, damaged cellular DNA (containing fragments and strand breaks) is separated from intact DNA, yielding a classic “comet tail” shape under the microscope. Extent of DNA damage is usually visually estimated by comet tail measurement; however, image analysis software is also available for measuring various parameters.

The OxiSelect™ Comet Assay is a fast and sensitive kit for the measurement of cellular DNA damage. Each kit provides sufficient reagents to perform up to 15 assays.

## **Assay Principle**

Cell Biolabs' OxiSelect™ Comet Assay is a single cell gel electrophoresis assay (SCGE) for simple evaluation of cellular DNA damage. First, individual cells are mixed with molten agarose before application to the OxiSelect™ Comet Slide. These embedded cells are then treated with a lysis buffer and alkaline solution, which relaxes and denatures the DNA. Finally, the samples are electrophoresed in a horizontal chamber to separate intact DNA from damaged fragments. Following electrophoresis, the samples are dried, stained with a DNA dye, and visualized by epifluorescence microscopy. Under these conditions, the damaged DNA (containing cleavage and strand breaks) will migrate further than intact DNA and produce a “comet tail” shape (see Figure 1). Each kit provides sufficient reagents to perform up to 15 assays.



**Figure 1: Comet Assay Principle**

## **Related Products**

1. STA-320: OxiSelect™ Oxidative DNA Damage ELISA Kit (8-OHdG Quantitation)
2. STA-321: OxiSelect™ DNA Double-Strand Break (DSB) Staining Kit
3. STA-324: OxiSelect™ Oxidative DNA Damage Quantitation Kit (AP sites)
4. STA-325: OxiSelect™ Oxidative RNA Damage ELISA Kit (8-OHG Quantitation)
5. STA-351: OxiSelect™ Comet Assay Kit (3-Well Slides)

## **Kit Components**

1. OxiSelect™ 3-Well Comet Slides (Part No. STA-352): Five slides.
2. OxiSelect™ Comet Agarose (Part No. 235002): One sterile 15 mL bottle.
3. Vista Green DNA Dye, 10000X (Part No. 235003): One 5 µL vial.
4. EDTA Solution, 500 mM (Part No. 235004): One 50 mL bottle.
5. 10X Lysis Solution (Part No. 235005): One 20 mL bottle.

## **Materials Not Supplied**

1. NaCl powder
2. NaOH pellets
3. 10 N NaOH for pH adjustment
4. DMSO (optional)
5. 70% Ethanol
6. TE Buffer (10 mM Tris, pH 7.5, 1 mM EDTA)
7. PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>)
8. EDTA (disodium salt)
9. DI H<sub>2</sub>O

## **Storage**

Upon receipt, store the Vista Green DNA Dye at -20°C. Store all other kit components at room temperature.

## **Preparation of Reagents**

- OxiSelect™ Comet Agarose: Heat the Comet Agarose bottle at 90-95°C in a water bath for 20 minutes, or until agarose liquefies. Transfer the bottle to a 37°C water bath for 20 minutes and maintain until needed.
- Vista Green DNA Dye: Prepare a 1X Vista Green DNA Staining Solution by diluting the provided stock 1:10000 in TE Buffer (10 mM Tris, pH 7.5, 1 mM EDTA). The solution can be stored at 4°C for up to 3 weeks, protected from light.
- Lysis Buffer: To prepare 100 mL of 1X Lysis Buffer

NaCl	14.6 g
EDTA Solution (provided)	20.0 mL
10X Lysis Solution (provided)	10.0 mL
DMSO	10.0 mL (optional for heme containing samples)
DI H <sub>2</sub> O	Adjust volume to 90 mL

Mix thoroughly to dissolve NaCl. Slowly adjust the Lysis Buffer to pH 10.0 with 10 N NaOH, then QS to 100 mL with DI H<sub>2</sub>O. Chill Lysis Buffer to 4°C before use.

*Note: Buffer will appear cloudy at room temperature, but will clear at 4°C. pH will also remain ~10.0.*

- Alkaline Solution: To prepare 100 mL of Alkaline Solution

NaOH	1.2 g
EDTA Solution (provided)	0.2 mL
DI H <sub>2</sub> O	Adjust volume to 100 mL

Mix thoroughly to dissolve NaOH. Chill Alkaline Solution to 4°C before use.

- Electrophoresis Running Solution: Choose the appropriate electrophoresis solution based on the desired running conditions and assay sensitivity. TBE is preferred for analysis of apoptosis and enables use of the tail length, rather than the tail moment, for data analysis. TBE electrophoresis will detect single-stranded and double-stranded DNA breaks, and may detect a few AP sites. Alkaline electrophoresis is more sensitive and will detect smaller amounts of DNA damage. Alkaline electrophoresis will detect single-stranded and double-stranded DNA breaks, the majority of AP sites, and alkali labile DNA adducts.

To prepare 1 L of Electrophoresis Solution

### 1. TBE Electrophoresis Solution

Tris Base	10.8 g
Boric Acid	5.5 g
EDTA (disodium salt)	0.93 g
DI H <sub>2</sub> O	Adjust volume to 1 L

Mix thoroughly to dissolve solids. Chill TBE Running Solution to 4°C before use.

### OR 2. Alkaline Electrophoresis Solution (300 mM NaOH, pH >13, 1 mM EDTA)

NaOH	12.0 g
EDTA Solution (provided)	2.0 mL
DI H <sub>2</sub> O	Adjust volume to 1 L

Mix thoroughly to dissolve NaOH. Chill Alkaline Running Solution to 4°C before use.

## **Special Precautions**

To avoid ultraviolet light damage to cell samples, perform the assay under low/dim light conditions.

## **Preparation of Samples and Slides**

1. Prepare Lysis Buffer, Alkaline Solution, and Electrophoresis Running Solution (see Preparation of Reagents) prior to performing the assay. Chill all solutions to 4°C thoroughly.
2. Heat OxiSelect™ Comet Agarose to 90-95°C in a water bath for 20 minutes, or until agarose liquefies. Cool the agarose by transferring the bottle to a 37°C water bath for 20 minutes.
3. Prepare cell samples, including controls, as follows:
  - Suspension Cells: Centrifuge cells at 700 x g for 2 minutes and discard supernatant. Wash cell pellet once with ice-cold PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), centrifuge, and discard the supernatant. Finally, resuspend the cells at 1 x 10<sup>5</sup> cells/mL in ice-cold PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>).
  - Adherent Cells: Gently remove cells from flask/dish by scraping with a rubber policeman. Transfer cell suspension to a conical tube and centrifuge at 700 x g for 2 minutes, discarding the supernatant. Wash cell pellet once with ice-cold PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), centrifuge, and discard the supernatant. Finally, resuspend the cells at 1 x 10<sup>5</sup> cells/mL in ice-cold PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>).
  - Tissue Preparation: Using dissection scissors, mince a small piece of tissue in 1-2 mL of ice-cold PBS containing 20 mM EDTA (without Mg<sup>2+</sup> and Ca<sup>2+</sup>). Allow the tissue/cell suspension to stand for 5 minutes before transferring the supernatant to a centrifuge tube; avoid transferring debris. Centrifuge, discarding the supernatant, and then resuspend the cells at 1 x 10<sup>5</sup> cells/mL in ice-cold PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>).
4. Combine cell samples with Comet Agarose (step 2) at 1:10 ratio (v/v), mix well by pipetting, and immediately transfer 75 µL/well onto the OxiSelect™ Comet Slide. **Ensure complete well coverage by spreading the solution over the well with the pipette tip.**

*Note: For multiple samples, warm the slide at 37°C, maintain suspensions at 37°C to avoid gelation. Titrate samples again just prior to slide addition*

5. Maintaining the slide horizontally, transfer the slide to 4°C in the dark for 15-30 minutes.
6. Carefully, transfer the slide to a small basin/container containing pre-chilled Lysis Buffer (~25 mL/slide). Immerse the slide in the buffer for 30-60 minutes at 4°C in the dark.
7. Carefully, aspirate the Lysis Buffer from the container and replace with pre-chilled Alkaline Solution (~25 mL/slide). Immerse the slide in the solution for 30 minutes at 4°C in the dark.

## **Assay Protocol**

### **I. TBE Electrophoresis**

1. Aspirate the Alkaline Solution from the container and replace with pre-chilled TBE Electrophoresis Solution. Immerse the slide for 5 minutes, and then repeat once more.


2. Maintaining the slide horizontally, carefully transfer the slide to a horizontal electrophoresis chamber. Fill the chamber with cold TBE Electrophoresis Solution until the buffer level covers the slide.
3. Apply voltage to the chamber for 10-15 minutes at 1 volt/cm (e.g., if the chamber electrodes are 35 cm apart, you would then apply 35 volts to the slide)
4. Maintaining the slide horizontally, carefully transfer the slide from the electrophoresis chamber to a clean, small basin/container containing pre-chilled DI H<sub>2</sub>O (~25 mL/slide). Immerse the slide for 2 minutes, aspirate, and then repeat twice more.
5. Aspirate the final water rinse and replace with cold 70% Ethanol for 5 minutes.
6. Maintaining the slide horizontally, remove the slide from the 70% Ethanol and allow slide to dry at 37°C for 30 min.
7. Once the agarose and slide are completely dry, add 100 µL/well of diluted Vista Green DNA Dye (see Preparation of Reagents). Incubate at room temperature for 15 minutes.
8. View slides by epifluorescence microscopy using a FITC filter.

## **II. Alkaline Electrophoresis**

1. Maintaining the slide horizontally, carefully transfer the slide from the Alkaline Solution to a horizontal electrophoresis chamber. Fill the chamber with cold Alkaline Electrophoresis Solution until the buffer level covers the slide.
2. Apply voltage to the chamber for 15-30 minutes at 1 volt/cm (e.g., if the chamber electrodes are 35 cm apart, you would then apply 35 volts to the slide). Additionally, adjust the volume of Alkaline Electrophoresis Solution to produce a current setting of 300 mA.
3. Maintaining the slide horizontally, carefully transfer the slide from the electrophoresis chamber to a clean, small basin/container containing pre-chilled DI H<sub>2</sub>O (~25 mL/slide). Immerse the slide for 2 minutes, aspirate, and then repeat twice more.
4. Aspirate the final water rinse and replace with cold 70% Ethanol for 5 minutes.
5. Maintaining the slide horizontally, remove the slide from the 70% Ethanol and allow slide to dry at 37°C for 30 min.
6. Once the agarose and slide are completely dry, add 100 µL/well of diluted Vista Green DNA Dye (see Preparation of Reagents). Incubate at room temperature for 15 minutes.
7. View slides by epifluorescence microscopy using a FITC filter.

## **Example of Results**


The following figures demonstrate typical OxiSelect™ Comet Assay Kit results. One should use the data below for reference only. This data should not be used to interpret actual results.



**Figure 2. Etoposide Treatment of Jurkat Cells.** Jurkat cells were untreated (left) or treated (right) with 20  $\mu$ M Etoposide for 4 hours before performing Comet Assay (alkaline electrophoresis conditions, 33 V/300 mA for 15 minutes).

### Calculation of Results

The DNA damage is quantified by measuring the displacement between the genetic material of the nucleus ('comet head') and the resulting 'tail'. Tail Moment and Tail DNA% are the two most common parameters to analyze Comet assay results. At least 50 -100 cells should be analyzed per sample. The Tail Moment has been suggested to be an appropriate index of induced DNA damage in considering both the migration of the genetic material as well as the relative amount of DNA in the tail.



**Figure 3: Typical Damaged DNA in Comet Assay.**

$$\text{Tail DNA\%} = 100 \times \frac{\text{Tail DNA Intensity}}{\text{Cell DNA Intensity}}$$

Tail Moment can be measured using one of the following methods:

(a) Olive Tail Moment = Tail DNA% x Tail Moment Length\*

(b) Extent Tail Moment = Tail DNA% x Length of Tail (see Figure 3)

A number of Comet Assay analysis software programs are commercially available, such as and Comet Assay IV (Perceptive Instruments) and CASPlab.

\*Tail Moment Length is measured from the center of the head to the center of the tail (see Figure 3)

### References

1. Ostling, O., and Johanson, K. J. (1984). *Biochem. Biophys. Res. Commun.* **123**, 291–298.
2. Singh, N. P., McCoy, M. T., Tice, R. R., and Schneider, E. L. (1988). *Exp. Cell. Res.* **175**, 184–191.
3. Olive, P. L., Banath, J. P., and Durand, R. E. (1990a). *Radiat. Res.* **122**, 86–94.



4. De Boeck, M., Touil, N., De Visscher, G., Vande, P. A., and Kirsch-Volders, M. (2000). *Mutat. Res.* **469**, 181–197.

### **Recent Product Citations**

1. Ni W. et al. (2023). An inducible long noncoding RNA, LncZFHX2, facilitates DNA repair to mediate osteoarthritis pathology. *Redox Biol.* doi: 10.1016/j.redox.2023.102858.
2. Tang, H. et al. (2023). Astaxanthin attenuated cigarette smoke extract-induced apoptosis via decreasing oxidative DNA damage in airway epithelium. *Biomed Pharmacother.* **167**:115471. doi: 10.1016/j.biopha.2023.115471.
3. Wong, H.T. et al. (2023). Inhibition of ATM-directed antiviral responses by HIV-1 Vif. *PLoS Pathog.* **19**(9):e1011634. doi: 10.1371/journal.ppat.1011634.
4. Hwang, Y.J. et al. (2023). StemRegenin 1 Mitigates Radiation-Mediated Hematopoietic Injury by Modulating Radioresponse of Hematopoietic Stem/Progenitor Cells. *Biomedicines.* **11**(3):824. doi: 10.3390/biomedicines11030824.
5. Elbadawi, M. et al. (2023). The Novel Artemisinin Dimer Isoniazide ELI-XXIII-98-2 Induces c-MYC Inhibition, DNA Damage, and Autophagy in Leukemia Cells. *Pharmaceutics.* **15**(4):1107. doi: 10.3390/pharmaceutics15041107.
6. Miu, B.A. et al. (2023). MRC-5 Human Lung Fibroblasts Alleviate the Genotoxic Effect of Fe-N Co-Doped Titanium Dioxide Nanoparticles through an OGG1/2-Dependent Reparatory Mechanism. *Int J Mol Sci.* **24**(7):6401. doi: 10.3390/ijms24076401.
7. Brabson, JP. Et al. (2023). Oxidized mC modulates synthetic lethality to PARP inhibitors for the treatment of leukemia. *Cell Rep.* **42**(1):112027. doi: 10.1016/j.celrep.2023.112027.
8. Calses, P.C. et al. (2023). TEAD Proteins Associate With DNA Repair Proteins to Facilitate Cellular Recovery From DNA Damage. *Mol Cell Proteomics.* **22**(2):100496. doi: 10.1016/j.mcpro.2023.100496.
9. Donne, R. et al. (2022). Replication stress triggered by nucleotide pool imbalance drives DNA damage and cGAS-STING pathway activation in NAFLD. *Dev Cell.* **57**(14):1728-1741.e6. doi: 10.1016/j.devcel.2022.06.003.
10. Zhu, Y. et al. (2022). Natural product preferentially targets redox and metabolic adaptations and aberrantly active STAT3 to inhibit breast tumor growth in vivo. *Cell Death Dis.* **13**(12):1022. doi: 10.1038/s41419-022-05477-2.
11. Yu, T. et al. (2022). Premature aging is associated with higher levels of 8-oxoguanine and increased DNA damage in the Polg mutator mouse. *Aging Cell.* **21**(9): e13669. doi: 10.1111/accel.13669.
12. Siemionow, M. et al. (2022). Long-Term Biodistribution and Safety of Human Dystrophin Expressing Chimeric Cell Therapy After Systemic-Intraosseous Administration to Duchenne Muscular Dystrophy Model. *Arch Immunol Ther Exp (Warsz).* **70**(1):20. doi: 10.1007/s00005-022-00656-7.
13. Hameed, F.H. & AL-Qadhi, H.I. (2022). Effect of Azithromycin on Sperm DNA of Male Rats. *Int. J. Drug Deliv. Technol.* **12**(2):594-597. doi: 10.25258/ijddt.12.2.22.
14. García-Cuellar, C.M. et al. (2022). Cetylpyridinium chloride inhibits human breast tumor cells growth in a no-selective way. *J Appl Biomater Funct Mater.* doi: 10.1177/22808000221092157.
15. Xiong, H. et al. (2022). IFN- $\gamma$  activates the tumor cell-intrinsic STING pathway through the induction of DNA damage and cytosolic dsDNA formation. *Oncoimmunology.* **11**(1):2044103. doi: 10.1080/2162402X.2022.2044103.

16. Albano, G.D. et al. (2022). Cadmium and Cadmium/BDE (47 or 209) Exposure Affect Mitochondrial Function, DNA Damage/Repair Mechanisms and Barrier Integrity in Airway Epithelial Cells. *Atmosphere*. **13**(2):201. doi: 10.3390/atmos13020201.
17. Fouché, T. et al. (2022). Ecotoxicological Effects of Aflatoxins on Earthworms under Different Temperature and Moisture Conditions. *Toxins*. **14**(2):75. doi: 10.3390/toxins14020075.
18. Zhang, L. et al. (2022). 53BP1 regulates heterochromatin through liquid phase separation. *Nat Commun*. **13**(1):360. doi: 10.1038/s41467-022-28019-y.
19. Chesnokova, V. et al. (2021). Local non-pituitary growth hormone is induced with aging and facilitates epithelial damage. *Cell Rep*. **37**(11):110068. doi: 10.1016/j.celrep.2021.110068.
20. Banerjee, D. et al. (2021). A non-canonical, interferon-independent signaling activity of cGAMP triggers DNA damage response signaling. *Nat Commun*. **12**(1):6207. doi: 10.1038/s41467-021-26240-9.
21. Yanuaryska, R.D. et al. (2021). Viability and DNA damage of buccal mucosa cells in patients exposed to panoramic X-ray. *Arch Orofac Sci*. **16**(Supp.1): 43–49. doi: 10.21315/aos2021.16.s1.8.
22. Talluri, S. et al. (2021). Dysregulated APOBEC3G causes DNA damage and promotes genomic instability in multiple myeloma. *Blood Cancer J*. **11**(10):166. doi: 10.1038/s41408-021-00554-9.
23. Wang, T. et al. (2021). The effects of glucose-6-phosphate dehydrogenase deficiency on benzene-induced hematotoxicity in mice. *Ecotoxicol Environ Saf*. **226**:112803. doi: 10.1016/j.ecoenv.2021.112803.
24. Ciminera, A.K. et al. (2021). Elevated glucose increases genomic instability by inhibiting nucleotide excision repair. *Life Sci Alliance*. **4**(10): e202101159. doi: 10.26508/lsa.202101159.
25. Hudita, A. et al. (2021). Bioinspired silk fibroin nano-delivery systems protect against 5-FU induced gastrointestinal mucositis in a mouse model and display antitumor effects on HT-29 colorectal cancer cells in vitro. *Nanotoxicology*. doi: 10.1080/17435390.2021.1943032.
26. Hung, S.Y. et al. (2021). Bavachinin Induces G2/M Cell Cycle Arrest and Apoptosis via the ATM/ATR Signaling Pathway in Human Small Cell Lung Cancer and Shows an Antitumor Effect in the Xenograft Model. *J Agric Food Chem*. doi: 10.1021/acs.jafc.1c01657.
27. Cho, K. et al. Suppressor of cytokine signaling 2 is induced in Huntington's disease and involved in autophagy. *Biochem Biophys Res Commun*. **559**:21-27. doi: 10.1016/j.bbrc.2021.04.089.
28. Cho, D.H. et al. (2021). Far-infrared irradiation inhibits breast cancer cell proliferation independently of DNA damage through increased nuclear Ca<sup>2+</sup>/calmodulin binding modulated-activation of checkpoint kinase 2. *J Photochem Photobiol B*. doi: 10.1016/j.jphotobiol.2021.112188.
29. Li, J. et al. (2021). Melatonin ameliorates cypermethrin-induced impairments by regulating oxidative stress, DNA damage and apoptosis in porcine Sertoli cells. *Theriogenology*. **167**:67-76. doi: 10.1016/j.theriogenology.2021.03.011.
30. Li, M.Z. et al. (2021). Discovery of MTR-106 as a highly potent G-quadruplex stabilizer for treating BRCA-deficient cancers. *Invest New Drugs*. doi: 10.1007/s10637-021-01096-4.

## **Warranty**

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

## **Contact Information**

Cell Biolabs, Inc.  
7758 Arjons Drive  
San Diego, CA 92126  
Worldwide: +1 858-271-6500  
USA Toll-Free: 1-888-CBL-0505  
E-mail: [tech@cellbiolabs.com](mailto:tech@cellbiolabs.com)  
[www.cellbiolabs.com](http://www.cellbiolabs.com)

©2008-2024: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.