Protocol for use No. 74: upcyte® HUVECs & Medium



(v1.1 - 11.2021)

When breaking the sealed bag of upcyte® cells you are explicitly accepting the terms of the limited use labellicense provided with the purchase of the cells. IT IS STRICTLY PROHIBITED TO EXPAND THE CELLS.

Unless indicated otherwise, upcyte technologies products and services are for research purpose only. Do notuse for diagnostic or therapeutic applications.

This protocol describes how to prepare upcyte® Endothelial Cell Culture Medium and how to culture upcyte® Human Umbilical Vein Endothelial Cells (HUVECs).

Required products for upcyte® HUVEC culture

(A) Endothelial Cell Culture Medium

Our Endothelial Cell Culture Medium is suitable for culture of primary human endothelial cells and upcyte® HUVECs. It consists of the following components:

- 500 mL basal medium
- 25 mL fetal bovine serum (FBS)
- 29.2 mL L-glutamine
- 2.5 mL supplement A
- 1.5 mL supplement B

Store the basal medium protected from light at 2-8 °C. Do not freeze the basal medium. Store FBS, Lglutamine and the supplements A/B at -20 °C. If stored properly, single components are

stable as indicated by the expiry date on the label.

For preparation of the complete medium, thaw FBS, L-glutamine and supplements A/B at room temperature. Transfer the entire content of each compound to the basal medium. Close the bottle, swirl gently and store at 2-8°C. After combining all components, the medium is stable for 6 weeks.

(B) upcyte® Human Umbilical Vein Endothelial Cells

Upon arrival, store upcyte® HUVECs in liquid or vapor phase nitrogen. They should not be stored at -70°C. After thawing, we recommend a preculture of 4-5 days to ensure sufficient recovery of

upcyte technologies GmbH Osterfeldstraße 12-14 22529 Hamburg, Germany

the cells.

(C) Additional products not supplied:

- PBS without Ca²⁺ or Mg²⁺
- Trypsin-EDTA (0.05% Trypsin / 0.02% EDTA)
- cell culture vessels and disposables

Culture protocol

Step 1: thawing of cells

- 1. Remove the cells from the storage tank.
- 2. Thaw cells at 37°C in a water bath for approximately 120 sec. A small piece of ice should still be visible. Disinfect the vial using 70% ethanol and transfer the vial to a laminar flow-hood.
- 3. Transfer the thawed cell suspension into a tube containing 10 mL of pre-warmed medium.
- 4. Wash the cryovial once with 1 mL of medium to maximize cell recovery.
- 5. Pellet the cells by centrifuging at 250 x q for 5 min at RT.
- 6. Aspirate the supernatant without disrupting the pellet.
- 7. Add an appropriate volume of pre-warmed medium to the pellet (2 mL per 2E6 cells) and carefully resuspend the cells.
- 8. Determine the cell number e.g. by using a Neubauer haemocytometer.
- 9. Seed the cells at 1E4 cells per cm². Use ~0.2 mL medium per cm² surface area.
- 10. Incubate at 37 °C and 5% CO₂.
- 11. Change the medium every 2-3 days.
- 12. Split the cells when reaching 70-80% confluence, usually achieved within 4-5 days.


Step 2: plating of cells

- Pre-warm culture medium to 37°C.
- Aspirate the medium from the culture flas

- 3. Wash the flask once with PBS.
- Aspirate PBS and add an appropriate volume of 0.05% trypsin/EDTA solution per culture dish (~20 μL per cm²).
- 5. Incubate 3-5 min at 37°C until most of the cells are rounded up (check under the microscope).
- 6. Gently tap the cell culture vessel to detach Remaining cells.
- 7. Add an equal volume of pre-warmed medium and rinse the culture surface.
- 8. Transfer the complete suspension to a tube and centrifuge for 5 min 250 x *g* at RT.
- 9. Discard supernatant and carefully resuspend the pellet in medium.
- 10. Determine the cell number as described in the previous section.
- 11. Seed 20 000 cells per cm² in your format of choice.

Morphology of upcyte® HUVECs

Cells initially show marked spreading but become smaller and compact at higher cell densities.



Morphology of upcyte® HUVECs one day after thawing


Unless indicated otherwise, upcyte technologies products and services are for research purpose only. Do not use for diagnostic or therapeutic applications.

Protocol for use No. 74: upcyte® HUVECs & Medium



When breaking the sealed bag of upcyte® cells you are explicitly accepting the terms of the limited use labellicense provided with the purchase of the cells. IT IS STRICTLY PROHIBITED TO EXPAND THE CELLS.

Unless indicated otherwise, upcyte technologies products and services are for research purpose only. Do notuse for diagnostic or therapeutic applications.



Morphology of upcyte® HUVECs after 3 days of culture

Product information

Product	Supplements/Components	Product number
Endothelial Cell Culture Medium	Basal Medium (500mL)	MEC003
	• FBS (25 mL)	
	L-Glutamine (29.2 mL)	
	Supplement A (2.5 mL)	
	Supplement B (1.5 mL)	
upcyte®	• 1 vial (2x10 ⁶ cells, cryopreserved)	CUV001
Human Umbilical Vein Endothelial		
Cells (HUVECs)		

Purchaser Notification

Limited Use Label License (upcyte technologies) This cellular product (upcyte® cells) generated using the upcyte® technology is provided under an intellectual property license from upcyte technologies GmbH. The transfer of the upcyte® cells is conditioned on the buyer using the purchased product solely in research conducted by the buyer, and the buyer must not use, sell or otherwise transfer this product, its components or any data generated with the product

- (a) for diagnostic, therapeutic or prophylactic purposes; or
- (b) for resale; or
- (c) for the production of therapeutic, diagnostic, prophylactic or any other products; or
- (d) to provide a service to deliver information or materials to a third party.

This intellectual property license is timely limited and granted to the buyer for 12 months after delivery of the upcyte® cells. After this time the buyer has to destroy remaining upcyte® cells and populations thereof or obtain a follow-up license. For information on purchasing a license to this upcyte technologies product for other purposes, for an extended time frame and/or for commercial purposes, contact upcyte technologies GmbH, Osterfeldstraße 12-14, 22529 Hamburg, Germany or info@upcyte.com.