

Taq 2X PCR

Master Mix

Catalog: RK20602

Size: 100 RXN / 500 RXN

Concentration: 2X

Components:

Taq 2X PCR Master Mix

RM20350

Product Description

Taq DNA Polymerase possesses a $5' \rightarrow 3'$ polymerase activity and 3' adenine (A) addition activity.

Taq 2X Master Mix is an optimized ready-to-use solution containing *Taq* DNA Polymerase, dNTPs, MgCl₂, KCI and stabilizers. It is ideally suited to routine PCR applications on various templates including pure DNA solutions, bacterial colonies, and cDNA products. It can amplify up to 4 kb from complex genomic DNA or up to 5 kb from lambda DNA. Applicable to PCR, colony PCR and primer extension.

Storage Temperature: -20 ℃

Heat Inactivation: No

5' - 3' Exonuclease: Yes

3' - 5' Exonuclease: No

Strand Displacement: +

Resulting Ends: Single-base 3 'Overhangs

Error Rate: ~ 285x10⁻⁶ bases

1X Master Mix Composition:

10 mM Tris-HCl, 1.5 mM MgCl $_2$, 50 mM KCl, 0.08% IPGAL 630, 0.05% Tween 20, pH8.6@25 ${\rm C}$; 200 μ M dNTPs, 5% Glycerol, 25 U/ml $\it Taq$ DNA Polymerase.

Instructions

Reaction setup:

We recommend assembling all reaction components on ice and quickly transferring the reactions to a thermocycler preheated to the denaturation temperature (95 $^{\circ}$ C).

Take 25 µl /50 µl system as an example.

Composition	25 µl	50 µl	Final Conc.
Nuclease-free H ₂ O	to 25 μl	to 50 μl	
10 μM Forward	0.5 μl	1 μl	0.2μΜ
Primer			$(0.05{\sim}1~\mu M)$
10 μM Reverse	0.5 μl	1 μl	0.2μΜ
Primer			$(0.05{\sim}1~\mu M)$
Template DNA	variable	variable	<1 μg/50 μl
Taq 2x PCR Master	12.5 µl	25 μΙ	1X
Mix			

Incubated in a thermocycler as the below program:

Temperature	Time	Cycles	
95 ℃	30s	1	
95 ℃	15-30s		
45-68 ℃	15-60s	30	
68 ℃	1 kb/min		
68 ℃	5min	1	
4-10 ℃	∞		

General Guidelines:

1. Template:

Use of high quality, purified DNA templates greatly enhances the success of PCR. Recommended amounts of DNA template for a 50 μ l reaction are as follows:

DNA	Amount
Genomic	1 ng-1 μg
Plasmid or viral	1 pg-1 ng

Order: order@abclonal.com

2. Primers:

Oligonucleotide primers are generally 20–40 nucleotides in length and ideally have a GC content of 40–60%. Computer programs such as Primer3 (http://frodo.wi.mit.edu/primer3) can be used to design or analyze primers. The final concentration of each primer in a reaction may be 0.05–1 μ M, typically 0.1–0.5 μ M.

3. Mg++ and additives:

Mg++ concentration of 1.5–2.0 mM is optimal for most PCR products generated with *Taq* DNA Polymerase. The final Mg++ concentration in 1X *Taq* PCR Master Mix is 1.5 mM. This supports satisfactory amplification of most amplicons. However, Mg++ can be further optimized in 0.5 or 1.0 mM increments using MgCl₂.

Amplification of some difficult targets, like GC-rich sequences, may be improved with additives, such as DMSO or formamide.

4. Denaturation:

An initial denaturation of 30 seconds at 95 $^{\circ}$ C is sufficient for most amplicons from pure DNA templates. For difficult templates such as GC-rich sequences, a longer denaturation of 2-4 minutes at 95 $^{\circ}$ C is recommended prior to PCR cycling to fully denature the template. With colony PCR, an initial 5 minutes denaturation at 95 $^{\circ}$ C is recommended.

During thermocycling a 15–30 second denaturation at 95 $^{\circ}$ C is recommended.

5. Annealing:

The annealing step is typically 15–60 seconds. Annealing temperature is based on the T_m of the primer pair and is typically 45–68 °C. Annealing temperatures can be optimized by doing a temperature gradient PCR starting 5 °C below the calculated T_m .

When primers with annealing temperatures above 65 $^{\circ}$ C are used, a 2-step PCR protocol is possible.

6. Extension:

The recommended extension temperature is $68 \, \mathbb{C}$. Extension times are generally 1 minute per kb. A final extension of 5 minutes at $68 \, \mathbb{C}$ is recommended.

7. Cycle number:

Generally, 25–35 cycles yields sufficient product. Up to 45 cycles may be required to detect low-copy-number targets.

8. 2-step PCR:

When primers with annealing temperatures above 65 $^{\circ}$ C are used, a 2-step thermocycling protocol is possible.

Thermocycling conditions for a routine 2-step PCR:

Temperature	Time	Cycles
95 ℃	30s	1
95 ℃ 65-68 ℃	15-30s 1kb/min	30
65-68 ℃ 4-10 ℃	5min ∞	1

9. PCR product:

The PCR products generated using *Taq* DNA Polymerase contain dA overhangs at the 3 ´-end; therefore the PCR products can be ligated to dT/dU-overhang vectors.