

Mouse anti-EBV-EBNA1, clone EBS-I-024 (Monoclonal)

Clone no. EBS-I-024

MONOSAN

Product name	Mouse anti-EBV-EBNA1, clone EBS-I-024 (Monoclonal)
Host	Mouse
Applications	ELISA, IF, IHC-P, WB
Species reactivity	EBV
Conjugate	-
Immunogen	EBNA1 synthetic peptide (AA-424-460)
Isotype	IgG1-K
Clonality	Monoclonal
Clone number	EBS-I-024
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-EBV-EBNA1, clone EBS-I-024 (Monoclonal)

Clone no. EBS-I-024

MONOSAN

Additional info

EBS-I-024 reacts with a conserved domain on EBNA1 (AA 430-442), present in all EBV isolates. EBNA1 is expressed in all viral latency stages in vitro and in vivo.

References

1. Chen, MR et al. J. Virol. 67(8): 4875-85 (1993)
2. Chen, MR et al. J. Biomed. Sci. 5(3): 173-9 (1998)
3. Sivachandra, N, et al, J. Virol. 86(1): 60-68 (2012)
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES