

Mouse anti-Vimentin, clone V9 (Monoclonal)

Clone no. V9

MONOSAN

---

Product name	Mouse anti-Vimentin, clone V9 (Monoclonal)
Host	Mouse
Applications	IHC-P (1:100), IHC-fr (1:100), WB, FC
Species reactivity	human
Conjugate	-
Immunogen	Vimentin isolated from cow lens
Isotype	IgG1
Clonality	Monoclonal
Clone number	V9
Size	1 ml
Concentration	25 ug/ml
Format	-
Storage buffer	Culture medium with 0.01M Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Vimentin, clone V9 (Monoclonal)

Clone no. V9

MONOSAN

**Additional info**

The antibody does stain all cells that contain vimentin. Only frozen section can be used. The antibody only stains the 57kDa vimentin band in immunoblots performed on a lysate of human lymphocytes.

**References**

1. Muijen GNP van (1984) Am J Pathol 116, 363-369
2. Debus E et al. (1984) Eur J Cell Biol 34, 137-143
3. Muijen G van (1986) Exp Cell Res 162; 97-113
4. Corver WE et al. (1994) Cytometry 15, 117-128
5. BroekJCM van den et al. (2000) Appl Immunohis & Mol Morphology 8, 316-321

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**