pMXs-Cas9-Puro Retroviral Vector

CATALOG NUMBER: RTV-5106 STORAGE: -80°C

QUANTITY AND CONCENTRATION: 100 µL of bacterial glycerol stock

Background

Retroviruses are efficient tools for delivering heritable genes into the genome of dividing cells. Cell Biolabs' pMXs-Puro retroviral vector is based on Moloney murine leukemia virus (MMLV). The vector provides the viral package signal, transcription and processing elements, and MCS for cloning of a target gene. The viral *env* gene, produced by the package cell line, encodes the envelope protein, which determines the viral infectivity range. Transfection into a package cell line produces high-titer, replication-incompetent viruses. In addition to transfer and expression of exogenous genes in mammalian cells, recently, retroviruses have been used to express silencing RNAs (siRNA) to decrease the expression of target genes both *in vitro* and *in vivo*.

The vector contains the ampicillin-resistance gene, MMLV LTRs, package signal, *S. pyogenes* Cas9 gene with NLSs (Nuclear Localization Signal), puromycin-resistance gene for the growth of infected mammalian cells to select stable cell lines (Figure 1).

Cas9 (CRISPR associated protein 9) is an RNA-guided DNA endonuclease. This enzyme associates with the CRISPR (Clustered Regularly Interspaced Short Palindromic Repeats) adaptive immunity system in various types of bacteria including *Streptococcus pyogenes*. Cas9 is able to unwind foreign DNA (such as plasmid DNA or invading bacteriophage DNA) and then checks for sites complementary to the 20 base pair spacer region of the guide RNA. If the DNA substrate is complementary to the guide RNA, Cas9 cuts up invading DNA.



Figure 1. Map of pMXs-Cas9-Puro Retroviral Vector.



Safety Consideration

Remember that you will be working with samples containing infectious virus. Follow the recommended NIH guidelines for all materials containing BSL-2 organisms. Always wear gloves, use filtered tips and work under a biosafety hood.

References

- 1. Kitamura T., et al., (2003) Exp. Hematol. 31, 1007-1014.
- 2. Heler, R; Samai, P; Modell, J. W.; Weiner, C; Goldberg, G. W.; Bikard, D; Marraffini, L. A. (2015). *Nature*. **519** :199–202.
- Jinek, M.; Chylinski, K.; Fonfara, I.; Hauer, M.; Doudna, J. A.; Charpentier, E. (2012). Science. 337: 816–821.
- 4. Mali, Prashant; Esvelt, Kevin M; Church, George M (2013). Nature Meth. 10: 957-963.
- 5. Mali, Prashant; Aach, John; Stranges, P Benjamin; Esvelt, Kevin M; Moosburner, Mark; Kosuri, Sriram; Yang, Luhan; Church, George M (2013). *Nature Biotech*. **31**: 833-838.
- 6. Gilbert, Luke A.; Larson, Matthew H.; Morsut, Leonardo; Liu, Zairan; Brar, Gloria A.; Torres, Sandra E.; Stern-Ginossar, Noam; Brandman, Onn; Whitehead, Evan H.; Doudna, Jennifer A.; Lim, Wendell A.; Weissman, Jonathan S.; Qi, Lei S. (2013). *Cell*. **154**: 442-451.
- 7. Zheng Q, Cai X, Tan MH, Schaffert S, Arnold CP, Gong X, Chen C-Z, and Huang S. (2014) BioTechniques **57**:115-124.

License Information

This product is licensed from the University of Tokyo.

Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

This product is for RESEARCH USE ONLY; not for use in diagnostic procedures.

Contact Information

Cell Biolabs, Inc. 7758 Arjons Drive San Diego, CA 92126 Worldwide: +1 858-271-6500 USA Toll-Free: 1-888-CBL-0505 E-mail: <u>tech@cellbiolabs.com</u> www.cellbiolabs.com

©2018: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

