

## Neo Glucose Solution (250 g/L)

#Cat: NB-58-0068

Size: 50ml

### General Information

D-(+)-Glucose, a common natural sugar, is involved in processes such as energy production, glycosylation and formation of glycans that provide structure to cells.

Glucose is used as a supplement for cell culture and in numerous cellular processes and molecular biology applications.

Applications:

- Additional component for cell culture media
- Applicable for batch, fed-batch & perfusion culture addition

### Product Specifications

<b>Appearance</b>	Clear solution
<b>Recommended concentration for media systems</b>	> 2 g/L
<b>Storage and shelf life</b>	Store at +2°C – +8°C protected from light. Once opened store at 4°C and use within 2-4 weeks.
<b>Shipping conditions</b>	Ambient

### Formulation

Components	Concentration mg/L
D-(+)-Glucose	250000.00

### Precautions and Disclaimer

This product is for research use only.

### Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email ([info@neo-biotech.com](mailto:info@neo-biotech.com)) or phone (+33 9 77 40 09 09).