

Mouse anti-Growth Hormone, clone GH-45 (Monoclonal)

Clone no. GH-45

MONOSAN

Product name	Mouse anti-Growth Hormone, clone GH-45 (Monoclonal)
Host	Mouse
Applications	IHC-P, ICC, ELISA
Species reactivity	Human
Conjugate	-
Immunogen	Human growth hormone.
Isotype	IgG1
Clonality	Monoclonal
Clone number	GH-45
Size	0.1 mg
Concentration	0.5 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON3061

MONOSAN[®]

Mouse anti-Growth Hormone, clone GH-45 (Monoclonal)

Clone no. GH-45

MONOSAN

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES