



## Monoclonal anti-GFP

**NB-01-0045**

**Volume:** 100  $\mu$ l

**Clones:** 1GFP-2A3

**Immunogen:** eGFP protein

**Murine Isotype:** IgG1, k

**Specificity:** GFP, GFP-fusion proteins

**Applications:**

- ELISA
- Western Blot
- Immunofluorescence
- Immunoprecipitation

**Recommended dilutions:** 1/500-1/5000

**Storage:**

Store at  $-20^{\circ}\text{C}$  (up to 3 years)

Avoid repeat freezing and thawing cycles. Does not contain any preservative