



MONOCLONAL ANTI-Gal 4 (DNA Binding Domain)

NB-01-0043

Volume: 100 µl

Format: Crude Ascites fluid

Immunogen: Gal 4 DBD fusion protein from aa 1 to aa147

Clones: 2GV-3

Murine Isotype: IgG1, k

Specificity: Gal 4 DBD containing proteins, epitope located in the second half of the DBD from aa 74 to 147

Applications:

- ELISA
- Immunoblotting
- Immunocytochemistry on transfected cells
- Immunoprecipitation
- Supershift experiments
- Transcription assays: see reference

Recommended dilutions: 1/1000-1/5000

Storage:

Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles. Does not contain any preservative

References:

White & al, EMBO J.(1992) vol.11, n°6: 2229-2240