

## Streptavidin-HRP Kit

Ready-to-use Streptavidin Conjugated Horse Radish Peroxidase Solution

Cat # NB-23-00026-1

Size: 500ml liquid, Ready-to-use

Cat # NB-23-00026-2

Size: 110ml liquid, Ready-to-use

Cat # NB-23-00026-3

Size: 18ml liquid, Ready-to-use

### Intended Use:

This enzyme conjugate solution can be used to label biotinylated antibodies.

### Recommended Protocol

- 1) Apply two drop (100ul) or enough volume to cover the tissue section completely
- 2) Incubate at room temperature for 10 minutes.
- 3) Resume standard immunohistochemistry staining procedure or find the procedures on our web site.

### Storage

Store at 2-8°C. Return the reagent back to 2-8°C after each use.

*For reference only*

*For Research Use Only. Not for Diagnostic or Therapeutic Use.*