Product datasheet MON9007-1



## Mouse anti-Renal Cell Carcinoma, clone RC38 (Monoclonal)

Clone no. RC38 MONOSAN

Product name Mouse anti-Renal Cell Carcinoma, clone RC38 (Monoclonal)

**Host** Mouse

Applications IHC-Fr (1:20)

Species reactivity human

Conjugate -

Immunogen Renal cell carcinoma homogenate

lsotype lgG1

**Clonality** Monoclonal

Clone number RC38

Size 1 ml

Concentration 15 ug/ml

Format -

Storage buffer Culture medium with 0.01M Sodium Azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON9007-1



Mouse anti-Renal Cell Carcinoma, clone RC38 (Monoclonal)

Clone no. RC38 MONOSAN

## Additional info

The antibody stains over 95% of primary renal cell carcinomas and 60% of metastases of renal cell carcinoma and glomerular visceral epithelium and proximal tubules in normal kidney (cytoplasmic). It does not stain epithelial cells of non-renal malignancies.

References

1. Oosterwijk E et al. (1986). Am J Pathol 123, 301-309

2 Dijk J. van (1988) J Nucl Med 29, 1087-1082

3. -

4.

5.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES