



Monoclonal anti-GCN5

NB-01-0044

Volume: 100 µl

Clones: 2GC-2C11

Immunogen: synthetic peptide coupled to ovalbumin: aa 118-132 (CKCNGWKNPKPPTAP) of human GCN5-Long, identical to mouse GCN5-Long in this region. Not present in human GCN5-Short .

Murine Isotype: IgG1, k

Specificity: human and mouse GCN5-L

Applications:

- ELISA
- Western Blot
- Immunofluorescence

Recommended dilutions: 1/500-1/5000

Storage:

Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles. Does not contain any preservative