

# NeoRoll PR22 - Plate Rotator

## NB-12-4016

NeoRoll PR22 was designed according to the feedback of users and their expectations about how a device should perform. As a result, this device is built to offer exceptional performance which is complemented by unprecedented ease-of-use. It's AC synchronous motor delivers all the power that you need reliably and without requiring frequent maintenance. Designed to maintain a small profile, NeoRoll PR22 wraps all its features and high-performance components in a small package to make the best use of limited bench space.



### Product Features



Brushless DC Motor For Quiet And Maintenance-free Long Life.

Variable Speed Setting  
From 10-80 RPM  
With Increments And  
Accuracy Of +/- 1 RPM.



Microprocessor Controlled Plate Rotator For Intelligent Performance With Last Run Memory Feature.

Variable Timer Setting  
From 1 To 99 Minutes  
With An Infinite Mode



Multiple Plate Options For Endless Flexibility. Supplied With One Chosen Disc Plate

User-programmable Feature Allows The User To Set Up To 9 User-defined Programs.





Fixed Speed Of  
30 RPM With An  
Accuracy Of +/- 2 RPM




Continuous/infinite Mode  
Of Operation For Keeping  
The Machine Running For  
Indefinite Lengths Of Time.

IP 21 Provides An  
Assurance Of  
Durability And Quality.



## Accessories



22 x 15 ml

4 x 50 ml + 12 x 15 ml

| Specifications                                     | NeoRoll PR22                       |
|--|------------------------------------|
| Capacity (supplied with any one chosen disc plate) | 12 x 15 ml + 4 x 50 ml, 22 x 15 ml |
| Motor Type   | AC Synchronous                     |
| Variable Speed                                     | Fix speed of 30 RPM                |
| Run Time   | Continuous / Infinite              |
| Speed Accuracy                                     | ± 2 RPM                            |
| Removable plate                                    | Yes                                |
| Noise Level  | Noiseless operation                |
| Weight (Capacity)                                  | 1 kg                               |
| Digital Display                                    | Yes                                |
| Dimension (WxDxH)                                  | 346 x 185 x 145 mm                 |
| Protection class                                   | IP 21                              |

## Ordering information

| Cat Number   | Details                                   |
|--------------|---|
| NB-12-4016-A | Plate rotator with 12 x 15 ml + 4 x 50 ml |
| NB-12-4016-B | Plate rotator with 15 ml clamp            |