

## **ENDO-Growth Media**

**Cat# NB-26-00280-500ML**

### **Description**

---

ENDO-Growth Media

### **Product Informations**

---

**Catalog #:** NB-26-00280-500ML

**Product Size:** 500ml

**Storage:** 2°C-8°C for up to 2 months  
-20°C long term storage

**Product Format:** Liquide

### **General Information**

---

ENDO-Growth Media contains two components freshly mixed for your convenience: 475 ml of ENDO-Basal media and 25 ml of ENDO-Growth supplement.

#### *ENDO-Basal Media*

Media containing central vitamins and amino acids, and 0.5% FBS.

#### *ENDO-Growth Supplement*

Contains the essential nutrients for endothelial cell proliferation; including 5% fetal bovine serum, recombinant growth factors and 1x penicillin and streptomycin

### **Shipping and Storage**

---

Product is shipped with ice packs. It can be stored at 2°C-8°C for up two months and for long term storage keep at -20°C.

## Protocol

---

1. Take the desired amount of ENDO-Growth Media need and warm to 37°C in a water bath or incubator
2. Avoid frequent temperature changes to the entire bottle of medium
3. The ENDO-Growth Media is stable for at least one month stored at 4°C
4. It is recommended that media is changed every day for regular cell culture

*For reference only*

***For Research Use Only. Not for Diagnostic or Therapeutic Use.***