

Personal Microplate Shaker - NeoShak Series

The Personal Microplate Shaker is an essential instrument for mixing individual micro plates of various kinds. The NeoShakers Series have a footprint as little as the plate dimension itself and provide mixing with orbital motion of 2 mm diameter and variable speed selection. There are 2 models one with variable speed only & the other (BL) with Variable Speed and Timer setting & USB connectivity.



NeoShak VT

Features :

- Tiny footprint - as small as the plate itself
- Universal fit - comes with interchangeable clamps for different micro plates
- Can be used inside fume hoods/cold rooms
- Comes with attachment for 20 x 1.5/2ml micro-tubes
- Optional USB connectivity for integration/remote operation
- Maintenance free BDLC Motor (NeoShak BL VT only)



NeoShak BL VT

Specifications :

Model	NeoShak VT	NeoShak BL VT
Cat#	NB-12-4005	NB-12-4006
Speed setting	200 to 1500 RPM	300 to 1800 RPM
Orbital diameter	2 mm	2 mm
Motor type	Stepper	Brushless DC motor
Timer setting	1 to 180min	Upto 99 mins + infinite mode
Display	---	Digital display of LED Speed & Time
Power supply	Universal (100 to 230V, 50/60 Hz)	
Dimensions (WxDxH)	127 x 85 x 86 mm	
Weight	1.4 Kg	
Standard attachments	<ul style="list-style-type: none"> • 20 x 1.5/2.0 ml Microtubes (Foam) • Replaceable clamps to use different types of standard micro plates 	



Foam Attachment



Tiny footprint



Turn Down Clamp

