



## Monoclonal anti-Estrogen receptors $\alpha$ and $\beta$

**NB-01-0027**

**Volume:** 100  $\mu$ l

**Clones:** 9ER $\beta$ -1F5

**Immunogen:** Recombinant human ER $\beta$  full-length protein

**Murine Isotype:** IgG2b, k

**Specificity:** Human and mouse estrogen receptors  $\alpha$  and  $\beta$ .

### **Applications:**

- ELISA
- Western Blot (natural and recombinant proteins)
- Immunofluorescence

**Recommended dilutions:** 1/500-1/5000

### **Storage:**

- Store at  $-20^{\circ}\text{C}$  (up to 3 years)
- Avoid repeat freezing and thawing cycles. Does not contain any preservative