A0076

Leader in Biomolecular Solutions for Life Science

# Integrin beta 3 (ITGB3/CD61) Rabbit pAb

Catalog No.: A0076



## **Basic Information**

**Observed MW** 97kDa

Calculated MW 87kDa

Category Polyclonal Antibody

Applications WB,ELISA

Cross-Reactivity Human,Mouse

### Background

The ITGB3 protein product is the integrin beta chain beta 3. Integrins are integral cellsurface proteins composed of an alpha chain and a beta chain. A given chain may combine with multiple partners resulting in different integrins. Integrin beta 3 is found along with the alpha IIb chain in platelets. Integrins are known to participate in cell adhesion as well as cellsurface mediated signalling.

# **Recommended Dilutions**

## Immunogen Information

WB

1:500 - 1:1000

**Gene ID** 3690 Swiss Prot P05106

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 609-788 of human ITGB3 (NP\_000203.2).

#### Synonyms

GT; GT2; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16; BDPLT24; Integrin beta 3 (ITGB3/CD61)

## Contact

# **Product Information**

www.abclonal.com

Purification Affinity purification

#### Storage

Source

Rabbit

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Isotype

lgG

## **Validation Data**

