



Monoclonal anti-ATF7 transcription factor

NB-01-0009

Clones: 1AT 2F10

Immunogen: recombinant human ATF7 protein (aa 1 to 494)

Murine Isotype: IgG1, k

Specificity: human ATF7 protein, epitope localized between aa 257 and 398.

Applications:

- ELISA
- Western Blot
- Immunofluorescence
- Immunoprecipitation
- gel supershift

Recommended dilutions: 1/500-1/5000

Storage:

- Store at -20°C (up to 3 years)
- Avoid repeat freezing and thawing cycles. Does not contain any preservative

References:

Hamard PJ, Dalbies-Tran R, Hauss C, Davidson I, Kedinger C, Chatton B. A functional interaction between ATF7 and TAF12 that is modulated by TAF4. *Oncogene*. 2005 May 12;24(21):3472-83.