



eIF4A1 antibody





Product Description

Mouse Monoclonal to eIF4A1.

Product Information

Code: NB-22-7036
Clonality Monoclonal

Host Mouse

Reactivity Human, Mouse, Rat

Applications WB, IF

ImmunogenSynthetic PeptideConjugationUnconjugated

Gene ID 1973 (Human); 13681 (Mouse);

Isotype lgG1

Dilution range WB 1:1000-3000; IF 1:100-200

Specificity The antibody detects endogenous eIF4A1 protein.

Swiss-Prot Key P60842 (Human); P60843 (Mouse)

Purification The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific

immunogen.

Form Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Concentration N/A Clone ID M8

Storage Store at -20°C. Avoid repeated freeze/thaw cycles.

Note For research use only.