


Bioss Easy-to-use IHC Kits

The data you can trust


Bioss easy-to-use IHC Kits are the perfect solution for quickly and conveniently obtaining high-quality IHC data. The kits contain all the reagents needed, not only including the primary and secondary antibodies, but also providing a positive control slide to confirm your staining. The optimized protocol is straightforward — simply use the reagents directly on your tissue sample and you are ready to gather your data for analysis.

What's in the box?

- | | | |
|---|--------------------------|--------------------------|
| ✓ Primary Antibody | ✓ Washing Buffer | ✓ Counter Stain |
| ✓ Antigen Retrieval Buffer | ✓ Secondary Antibody | ✓ Differentiation Buffer |
| ✓ Blocking Buffer | ✓ DAB Chromogen | ✓ Mounting Media |
| ✓ Endogenous Peroxidase Blocking Buffer | ✓ Positive Control Slide | ✓ Manual |



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using IHC0102M (Histone H3 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat cerebellum tissue slide using IHC0110R (S100B IHC Kit).

Cat No.	Protein Target	Primary Antibody Clonality
IHC0101H	Human GFAP	Rabbit Polyclonal
IHC0101M	Mouse GFAP	Rabbit Polyclonal
IHC0101R	Rat GFAP	Rabbit Polyclonal
IHC0102H	Human Histone H3	Rabbit Polyclonal
IHC0102M	Mouse Histone H3	Rabbit Polyclonal
IHC0102R	Rat Histone H3	Rabbit Polyclonal
IHC0105H	Human MYBL2	Rabbit Polyclonal
IHC0105M	Mouse MYBL2	Rabbit Polyclonal
IHC0105R	Rat MYBL2	Rabbit Polyclonal
IHC0106H	Human Calbindin	Rabbit Polyclonal
IHC0106M	Mouse Calbindin	Rabbit Polyclonal
IHC0106R	Rat Calbindin	Rabbit Polyclonal
IHC0107H	Human MBP	Rabbit Polyclonal
IHC0107M	Mouse MBP	Rabbit Polyclonal
IHC0107R	Rat MBP	Rabbit Polyclonal
IHC0108H	Human Peripherin	Rabbit Polyclonal
IHC0108M	Mouse Peripherin	Rabbit Polyclonal
IHC0108R	Rat Peripherin	Rabbit Polyclonal
IHC0110H	Human S100B	Mouse Monoclonal
IHC0110M	Mouse S100B	Mouse Monoclonal
IHC0110R	Rat S100B	Mouse Monoclonal
IHC0109H	Human GLP-I	Mouse Monoclonal
IHC0109M	Mouse GLP-I	Mouse Monoclonal
IHC0109R	Rat GLP-I	Mouse Monoclonal
IHC0104H	Human Glucagon	Rabbit Polyclonal
IHC0104M	Mouse Glucagon	Rabbit Polyclonal
IHC0104R	Rat Glucagon	Rabbit Polyclonal
IHC0114H	Human Alpha Smooth Muscle Actin	Mouse Monoclonal
IHC0114M	Mouse Alpha Smooth Muscle Actin	Mouse Monoclonal
IHC0114R	Rat Alpha Smooth Muscle Actin	Mouse Monoclonal
IHC0113H	Human Vimentin	Mouse Monoclonal
IHC0103M	Mouse TFF3	Rabbit Polyclonal
IHC0103R	Rat TFF3	Rabbit Polyclonal
IHC0112H	Human Cytokeratin 7	Mouse Monoclonal

CliniSciences Group

Austria



Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>

Belgium



Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>

Denmark



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>

Finland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>

France



Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>

Germany



Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>

Iceland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>

Ireland



Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>

Italy



Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>

Netherlands



Company: CliniSciences B.V.
Address: Kraaijenhoffstraat 137A,
1018RG Amsterdam, - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>

Norway



Company: CliniSciences AS
Address: c/o MerVerdi Munkerudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>

Poland



Company: CliniSciences sp.z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>

Portugal



Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>

Spain



Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 916 750 700
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>

Sweden



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>

Switzerland



Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>

UK



Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>

USA



Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>