

Crafted Exclusively for Pharma

IS Series Cytokine/ Biomarker ELISA Kits

CUSABIO's IS (Industry Standard) series of products, designed for industrial customers such as pharmaceutical companies, include cytokine and biomarker detection ELISA kits. Leveraging extensive R&D experience and independent development platforms, we screen for high-performance raw and auxiliary materials and adhere to strict production and quality control standards to ensure product excellence.

+ Application Areas

- CRS detection
- Residue detection in semi-finished drug processes
- Cell function evaluation in pharmacokinetic experiments
- Development quality control and pharmacodynamic testing of cell therapy drugs



+ Product Advantages



High-performance raw materials are strictly screened to provide ready-to-use finished products based on the requirements of industrial customers.



Inventory could be flexibly increased based on customized needs of pharmaceutical companies.



8 quality control parameters, strict quality control performed by professional teams.



Refined production quality control to minimize variation.

Human Interleukin 10,IL-10 ELISA KIT as an example:

Product Types	+ IS series industrial special products CSB-E04593h-IS						Conventional scientific research products CSB-E04593h					
Sample Types	serum, plasma, tissue homogenates, cell culture supernates						serum,plasma, tissue homogenates					
Detection Range	3.9 pg/mL-250 pg/mL						31.25 pg/mL- 2000 pg/mL					
Sensitivity	0.975 pg/mL						7.8 pg/mL					
Standard Curve	Standard Pg/mL 2.9010 250 2.1421 125 1.5607 62.5 1.0693 31.2 0.6319 15.6 0.3711 7.8 0.2450 3.9 0.1462 0						Standard Pg/mL 2.5902 2000 2.0821 1000 1.4685 500 0.9354 250 0.5238 125 0.3046 62.5 0.2114 31.25 0.1054 0					
Precision				Intra-Assay Precision			Inter-Assay Precision					
	Sample	1	2	3	1	2	3	Sample	1	2	3	
	n	20	20	20	20	20	20	n	20	20	20	
	Mean(pg/mL)	35.83	33.975	36.452	33.565	36.973	36.037	Mean(pg/mL)	251.241	254.176	249.776	246.854 255.646 251.974
	SD	0.031	0.024	0.029	0.041	0.034	0.036	SD	0.028	0.024	0.027	0.04 0.037 0.035
	CV(%)	5.688	4.554	5.263	7.839	6.115	6.581	CV(%)	5.283	4.494	5.114	7.634 6.903 6.591

© Cytokines

BDNF	CCL18	CCL2	CCL3	CCL4	CCL5	CXCL10	CXCL5
HGF	hs-CRP	IFNG	IFN-γ	IL10	IL12	IL13	IL17A
IL17F	IL1B	IL1β	IL2	IL21	IL22	IL29	IL4
IL6	IL7	TGF-β1	TNF-α				

© Biomarker

CD81	ERBB2	gp130	GPNMB	ICAM1	INS	Kim-1	MMP3
PRF1	RAGE/AGER	SIRT1	THPO	TNFSF10	TNFsR- I	Tp53	

*Scan to view all products.

Distributed by:

CliniSciences Group

CliniSciences Group

Austria



Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>

Belgium



Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>

Denmark



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>

Finland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>

France



Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>

Germany



Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>

Iceland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>

Ireland



Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>

Italy



Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>

Netherlands



Company: CliniSciences B.V.
Address: Kraijenhoffstraat 137A,
1018RG Amsterdam, - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>

Norway



Company: CliniSciences AS
Address: c/o MerVerdi Munkerudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>

Poland



Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>

Portugal



Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>

Spain



Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 916 750 700
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>

Sweden



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>

Switzerland



Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>

UK



Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>

USA



Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>